

# Mobile EasyBed

Good for Long-term Care, Nursing Homes,  
and Home Health Care

Designed for Persons  
who are Elderly or with Disabilities



Home Care  Nursing Home



YM2000B Standard Model



Invention Patents : USA/Taiwan  
Utility Model Patents : Japan/China  
U.S. FDA Approved

Made in Taiwan

Mobile EasyBed :  
Electric Bed + Wheelchair +  
Transfer Device

## ★Product Highlights★



Relieve the burden on  
caregivers with  
height-Adjustable function.



The bed chair can lay  
down or sit up.



The side door can be  
opened which allows you  
to pull out the bed chair.



The bed chair can  
work in limited space  
and easily operate  
both in and outdoors.

## PRODUCT

### YM2000A Alternating Air Model



### YM2000B Standard Model



Mobile EasyBed can relieve the burden on caregivers, which reduces the risk of further injuries for both caregivers and patients.

## Mobile EasyBed Specification

Bed (Material: S10C, SPHC, Iron)	Bed Load	173 Kg
	Weight	170(±1) kg (for YM2000A) 167(±1) kg (for YM2000B)
	Dimension (LxWxH)	2,280 x 1,030 x 910~1,060 mm
	Height/Low (Sleeping Side)	740(±10) mm / 590(±10) mm
	Mattress (LxWxH)	1900 x 500 x 127 mm (for YM2000A) 1940 x 520 x 100 mm (for YM2000B)
	Caster Wheel Diameter	5"(125 mm) x 4PCs [With Brake for 4 Bed Legs] 3"(75 mm) x 8PCs [Without Brake · Active Frame]
BedChair (Wheel Chair)	Cushion (LxWxH)	460 x 500 x 127 mm (for YM2000A) 460 x 520 x 100 mm (for YM2000B)
	Allow chair back rising Up/Down Angles	0~72 degrees(±2)
	Allow chair foot rising Up/Down Angles	+16 degrees~-75 degrees
	Caster Wheel Diameter	5"(125 mm) x 4PCs [With Brake]
Control Box	Input Voltage	100~125V AC / 220~230V AC
	Output Voltage	24V DC
	Bed Frame Line Controller	2 keys for Bed UP/ Down with Voice Notice
	Bed chair Line Controller	6 Keys for Foot/ Back Up/ Down
	Water Proof Level	IP53
	Certification	IEC60601-1-1 / IEC60601-1-2
Portable Battery	Rechargeable Battery Output 24V Voltage	2x [1.3A/12V] Green LED Indicator while Activating

## Order Information

Model	Description		
YM2000AL		Bed Chair Moves out from <b>Left Side</b>	Alternating Air Mattress only for Bed Chair, with 18 Fully Alternating Independent Cells
YM2000AR		Bed Chair Moves out from <b>Right Side</b>	
YM2000BL		Bed Chair Moves out from <b>Left Side</b>	Standard Mattress (Jersey Fabric Material) High Density Sponge with Polymer Cotton
YM2000BR		Bed Chair Moves out from <b>Right Side</b>	
AM011		Alternating Air Mattress	Alternating Air Mattress for YM2000 Bed Chair Use Only
AM021		Standard Mattress	Standard Mattress for YM2000 Bed Chair Use Only Jersey Fabric Material with High Density Sponge



YOFA Biotechnology Co., Ltd.

<https://en.yofa-tech.com/>  
[MAIL:sales@yofa-tech.com](mailto:sales@yofa-tech.com)

Official Website



Product Video

